

Recombinant Human GSPT2

Catalog No: #GP13437



Package Size: #GP13437-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GSPT2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 429-628 amino acids of human GSPT2
Target Name	G1 to S phase transition 2
Other Names	GST2; ERF3B
Accession No.	Swissprot:Q8IYD1Gene Accession:BC036077
Uniprot	Q8IYD1
GeneID	23708;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a GTPase that belongs to the GTP-binding elongation factor family. The encoded protein is a polypeptide release factor that complexes with eukaryotic peptide chain release factor 1 to mediate translation termination. This protein may also be involved in mRNA stability.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only