

# Recombinant Human SHC4

Catalog No: #GP13438



Package Size: #GP13438-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SHC4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 431-630 amino acids of human SHC4
Target Name	SHC adaptor protein 4
Other Names	RaLP; SHCD
Accession No.	Swissprot:Q6S5L8Gene Accession:BC033907
Uniprot	Q6S5L8
GeneID	399694;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Activates both Ras-dependent and Ras-independent migratory pathways in melanomas. Contributes to the early phases of agrin-induced tyrosine phosphorylation of CHRN1.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only