

Recombinant Human FARP2

Catalog No: #GP13441



Package Size: #GP13441-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human FARP2 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 1-200 amino acids of human FARP2 |
| Target Name | FERM, ARH/RhoGEF and pleckstrin domain protein 2 |
| Other Names | FIR; FRG; PLEKHC3 |
| Accession No. | Swissprot:O94887Gene Accession:BC021301 |
| Uniprot | O94887 |
| GeneID | 9855; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

Functions as guanine nucleotide exchange factor that activates RAC1. May have relatively low activity. Plays a role in the response to class 3 semaphorins and remodeling of the actin cytoskeleton. Plays a role in TNFSF11-mediated osteoclast differentiation, especially in podosome rearrangement and reorganization of the actin cytoskeleton. Regulates the activation of ITGB3, integrin signaling and cell adhesion (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only