

Recombinant Human TRO

Catalog No: #GP13449



Package Size: #GP13449-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TRO
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human TRO
Target Name	trophinin
Other Names	MAGED3; MAGE-d3
Accession No.	Swissprot:Q12816Gene Accession:BC026914
Uniprot	Q12816
GeneID	7216;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a membrane protein that mediates cell adhesion between trophoblastic cells and the epithelial cells of the endometrium. The encoded protein participates in cell signalling during embryo implantation, and may also be involved in cancer formation. This gene is located near several other closely related genes on chromosome X. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only