

# Recombinant Human TSGA10

Catalog No: #GP13452



Package Size: #GP13452-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TSGA10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 499-698 amino acids of human TSGA10
Target Name	testis specific 10
Other Names	CT79; CEP4L; SPGF26
Accession No.	Swissprot:Q9BZW7Gene Accession:BC028366
Uniprot	Q9BZW7
GeneID	80705;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May play a role in the sperm tail fibrous sheath, a major sperm tail structure.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only