

# Recombinant Human KANSL1L

Catalog No: #GP13453



Package Size: #GP13453-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human KANSL1L
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 477-676 amino acids of human KANSL1L
Target Name	KAT8 regulatory NSL complex subunit 1 like
Other Names	MSL1v2; C2orf67
Accession No.	Swissprot:A0AUZ9Gene Accession:BC032837
Uniprot	A0AUZ9
GeneID	151050;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only