

Recombinant Human DPY19L3

Catalog No: #GP13456



Package Size: #GP13456-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DPY19L3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-102 amino acids of human DPY19L3
Target Name	dpy-19 like C-mannosyltransferase 3
Accession No.	Swissprot:Q6ZPD9Gene Accession:BC029162
Uniprot	Q6ZPD9
GeneID	147991;
Storage	-20~-80°C, pH 7.6 PBS

Background

Probable C-mannosyltransferase that mediates C-mannosylation of tryptophan residues on target proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only