

WNT8A Antibody

Catalog No: #43193

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WNT8A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total WNT8A protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human WNT8A
Target Name	WNT8A
Other Names	WNT8D
Accession No.	Swiss-Prot#: Q9H1J5 Gene ID: 7478
Uniprot	Q9H1J5
GeneID	7478;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43193 at dilution 1/20,

Background

The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family, and may be implicated in development of early embryos as well as germ cell tumors. Multiple alternatively spliced transcript variants have been found for this gene.

Note: This product is for in vitro research use only