

# Recombinant Human MIA

Catalog No: #GP13480



Package Size: #GP13480-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MIA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 25-131 amino acids of human MIA
Target Name	melanoma inhibitory activity
Other Names	CD-RAP
Accession No.	Swissprot:Q16674Gene Accession:BC005910
Uniprot	Q16674
GeneID	8190;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only