

Recombinant Human TPD52L3

Catalog No: #GP13481



Package Size: #GP13481-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TPD52L3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-132 amino acids of human TPD52L3
Target Name	tumor protein D52 like 3
Other Names	D55; hD55; NYDSP25
Accession No.	Swissprot:Q96J77Gene Accession:BC017589
Uniprot	Q96J77
GeneID	89882;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the tumor protein D52-like family of proteins. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. The encoded protein may play a role in spermatogenesis. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only