

Recombinant Human EIF1AX

Catalog No: #GP13486



Package Size: #GP13486-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EIF1AX
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-144 amino acids of human EIF1AX
Target Name	eukaryotic translation initiation factor 1A X-linked
Other Names	EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1
Accession No.	Swissprot:P47813Gene Accession:BC000793
Uniprot	P47813
GeneID	1964;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA_i and eIF3) to the 5' end of capped RNA.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only