

Recombinant Human CPLX3

Catalog No: #GP13496



Package Size: #GP13496-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CPLX3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-158 amino acids of human CPLX3
Target Name	complexin 3
Other Names	CPXIII; CPX-III; Nbla11589
Accession No.	Swissprot:Q8WVH0Gene Accession:BC018026
Uniprot	Q8WVH0
GeneID	594855;
Storage	-20~-80°C, pH 7.6 PBS

Background

Positively regulates a late step in synaptic vesicle exocytosis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only