


Description

Product Name	BAZ1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BAZ1B protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human BAZ1B
Target Name	BAZ1B
Other Names	WSTF; WBSCR9; WBSCR10
Accession No.	Swiss-Prot#: Q9UIG0 Gene ID: 9031
Uniprot	Q9UIG0
GeneID	9031;
Concentration	2.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:45-1:300

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #43197 at dilution 1/35,

Background

This gene encodes a member of the bromodomain protein family. The bromodomain is a structural motif characteristic of proteins involved in chromatin-dependent regulation of transcription. This gene is deleted in Williams-Beuren syndrome, a developmental disorder caused by deletion of multiple genes at 7q11.23.?

Note: This product is for in vitro research use only