

Recombinant Human ZG16

Catalog No: #GP13500



Package Size: #GP13500-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ZG16
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 17-167 amino acids of human ZG16
Target Name	zymogen granule protein 16
Other Names	JCLN; JCLN1; ZG16A
Accession No.	Swissprot:O60844Gene Accession:BC029149
Uniprot	O60844
GeneID	653808;
Storage	-20~-80°C, pH 7.6 PBS

Background

May play a role in protein trafficking. May act as a linker molecule between the submembranous matrix on the luminal side of zymogen granule membrane (ZGM) and aggregated secretory proteins during granule formation in the TGN.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only