

Recombinant Human MYDGF

Catalog No: #GP13507



Package Size: #GP13507-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MYDGF
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 32-173 amino acids of human MYDGF
Target Name	myeloid derived growth factor
Other Names	IL25; IL27; SF20; IL27w; C19orf10; R33729_1; EUROIMAGE1875335
Accession No.	Swissprot:Q969H8Gene Accession:BC010129
Uniprot	Q969H8
GeneID	56005;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only