

# Recombinant Human RNF208

Catalog No: #GP13512



Package Size: #GP13512-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RNF208
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 82-261 amino acids of human RNF208
Target Name	ring finger protein 208
Accession No.	Swissprot:Q9H0X6Gene Accession:BC016958
Uniprot	Q9H0X6
GeneID	727800;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only