

Recombinant Human RPL18

Catalog No: #GP13515



Package Size: #GP13515-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RPL18
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-188 amino acids of human RPL18
Target Name	ribosomal protein L18
Other Names	L18
Accession No.	Swissprot:Q07020Gene Accession:BC000374
Uniprot	Q07020
GeneID	6141;
Storage	-20~-80°C, pH 7.6 PBS

Background

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L18E family of ribosomal proteins that is a component of the 60S subunit. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only