

# Recombinant Human MARCKSL1

Catalog No: #GP13520



Package Size: #GP13520-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human MARCKSL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-195 amino acids of human MARCKSL1
Target Name	MARCKS like 1
Other Names	F52; MLP; MRP; MLP1; MACMARCKS
Accession No.	Swissprot:P49006Gene Accession:BC007904
Uniprot	P49006
GeneID	65108;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the myristoylated alanine-rich C-kinase substrate (MARCKS) family. Members of this family play a role in cytoskeletal regulation, protein kinase C signaling and calmodulin signaling. The encoded protein affects the formation of adherens junction. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on the long arm of chromosomes 6 and 10.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only