

Recombinant Human COMMD9

Catalog No: #GP13522



Package Size: #GP13522-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COMMD9
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-198 amino acids of human COMMD9
Target Name	COMM domain containing 9
Other Names	HSPC166; C11orf55; LINC00610
Accession No.	Swissprot:Q9P000Gene Accession:BC010892
Uniprot	Q9P000
GeneID	29099;
Storage	-20~-80°C, pH 7.6 PBS

Background

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes (PubMed:21778237). May down-regulate activation of NF-kappa-B (PubMed:15799966). Modulates Na⁺ transport in epithelial cells by regulation of apical cell surface expression of amiloride-sensitive sodium channel (ENaC) subunits (PubMed:23637203).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only