

Recombinant Human JOSD1

Catalog No: #GP13528



Package Size: #GP13528-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human JOSD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-202 amino acids of human JOSD1
Target Name	Josephin domain containing 1
Other Names	dJ508115.2
Accession No.	Swissprot:Q15040Gene Accession:BC015026
Uniprot	Q15040
GeneID	9929;
Storage	-20~-80°C, pH 7.6 PBS

Background

Deubiquitinates monoubiquitinated probes (in vitro). When ubiquitinated, cleaves 'Lys-63'-linked and 'Lys-48'-linked poly-ubiquitin chains (in vitro), hence may act as a deubiquitinating enzyme. May increase macropinocytosis and suppress clathrin- and caveolae-mediated endocytosis. May enhance membrane dynamics and cell motility independently of its catalytic activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only