

Recombinant Human MYEF2

Catalog No: #GP13534



Package Size: #GP13534-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MYEF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 366-576 amino acids of human MYEF2
Target Name	myelin expression factor 2
Other Names	MEF-2; MST156; myEF-2; MSTP156; HsT18564
Accession No.	Swissprot:Q9P2K5Gene Accession:BC014533
Uniprot	Q9P2K5
GeneID	50804;
Storage	-20~-80°C, pH 7.6 PBS

Background

Transcriptional repressor of the myelin basic protein gene (MBP). Binds to the proximal MB1 element 5'-TTGTCC-3' of the MBP promoter. Its binding to MB1 and function are inhibited by PURA (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only