

Recombinant Human PPIC

Catalog No: #GP13536



Package Size: #GP13536-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PPIC
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 30-212 amino acids of human PPIC
Target Name	peptidylprolyl isomerase C
Other Names	CYPC
Accession No.	Swissprot:P45877Gene Accession:BC002678
Uniprot	P45877
GeneID	5480;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. Similar to other PPIases, this protein can bind immunosuppressant cyclosporin A.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only