

Recombinant Human PRPF40A

Catalog No: #GP13539



Package Size: #GP13539-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PRPF40A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 743-957 amino acids of human PRPF40A
Target Name	pre-mRNA processing factor 40 homolog A
Other Names	HYP A; FBP11; FLAF1; FNBP3; HIP10; Prp40; FBP-11; HIP-10; NY-REN-6
Accession No.	Swissprot:O75400Gene Accession:BC027178
Uniprot	O75400
GeneID	55660;
Storage	-20~-80°C, pH 7.6 PBS

Background

Binds to WASL/N-WASP and suppresses its translocation from the nucleus to the cytoplasm, thereby inhibiting its cytoplasmic function (By similarity). Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape and migration. May play a role in cytokinesis. May be involved in pre-mRNA splicing.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only