

Recombinant Human GINS3

Catalog No: #GP13540



Package Size: #GP13540-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GINS3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-216 amino acids of human GINS3
Target Name	GINS complex subunit 3
Other Names	PSF3
Accession No.	Swissprot:Q9BRX5Gene Accession:BC005879
Uniprot	Q9BRX5
GeneID	64785;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein subunit of the GINS heterotetrameric complex, which is essential for the initiation of DNA replication and replisome progression in eukaryotes. Alternatively spliced transcript variants encoding distinct isoforms have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only