

## Recombinant Human SLC30A6

Catalog No: #GP13541



Package Size: #GP13541-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SLC30A6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 246-461 amino acids of human SLC30A6
Target Name	solute carrier family 30 member 6
Other Names	ZNT6; MST103; MSTP103
Accession No.	Swissprot:Q6NXT4Gene Accession:BC032525
Uniprot	Q6NXT4
GeneID	55676;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of a family of proteins that function as zinc transporters. This protein can regulate subcellular levels of zinc in the Golgi and vesicles. Expression of this gene is altered in the Alzheimer's disease brain plaques.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only