

Recombinant Human FAAP100

Catalog No: #GP13551



Package Size: #GP13551-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FAAP100
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 649-881 amino acids of human FAAP100
Target Name	FA core complex associated protein 100
Other Names	C17orf70
Accession No.	Swissprot:Q0VG06Gene Accession:BC021968
Uniprot	Q0VG06
GeneID	80233;
Storage	-20~-80°C, pH 7.6 PBS

Background

FAAP100 is a component of the Fanconi anemia (FA; MIM 277650) core complex and is required for core complex stability and FANCD2 (see MIM 227646) monoubiquitination (Ling et al., 2007 [PubMed 17396147]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only