

Recombinant Human PREPL

Catalog No: #GP13555



Package Size: #GP13555-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PREPL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 512-727 amino acids of human PREPL
Target Name	prolyl endopeptidase like
Other Names	CMS22
Accession No.	Swissprot:Q4J6C6Gene Accession:BC013193
Uniprot	Q4J6C6
GeneID	9581;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only