

# Recombinant Human RTP4

Catalog No: #GP13558



Package Size: #GP13558-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RTP4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-168 amino acids of human RTP4
Target Name	receptor transporter protein 4
Other Names	IFRG28; Z3CXXC4
Accession No.	Swissprot:Q96DX8Gene Accession:BC013161
Uniprot	Q96DX8
GeneID	64108;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Probable chaperone protein which facilitates trafficking and functional cell surface expression of some G-protein coupled receptors (GPCRs). Promotes functional expression of the bitter taste receptor TAS2R16 (PubMed:16720576). Also promotes functional expression of the opioid receptor heterodimer OPRD1-OPRM1 (By similarity).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only