

Recombinant Human PLEKHF2

Catalog No: #GP13563

Package Size: #GP13563-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PLEKHF2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-249 amino acids of human PLEKHF2
Target Name	pleckstrin homology and FYVE domain containing 2
Other Names	EAPF; PHAFIN2; ZFYVE18
Accession No.	Swissprot:Q9H8W4Gene Accession:BC011806
Uniprot	Q9H8W4
GeneID	79666;
Storage	-20~-80°C, pH 7.6 PBS

Background

May play a role in early endosome fusion upstream of RAB5, hence regulating receptor trafficking and fluid-phase transport. Enhances cellular sensitivity to TNF-induced apoptosis (PubMed:18288467).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only