

Recombinant Human SUSD3

Catalog No: #GP13570



Package Size: #GP13570-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SUSD3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-255 amino acids of human SUSD3
Target Name	sushi domain containing 3
Accession No.	Swissprot:Q96L08Gene Accession:BC014601
Uniprot	Q96L08
GeneID	203328;
Storage	-20~-80°C, pH 7.6 PBS

Background

May play a role in breast tumorigenesis by promoting estrogen-dependent cell proliferation, cell-cell interactions and migration.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only