

Recombinant Human ELF5

Catalog No: #GP13571



Package Size: #GP13571-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human ELF5 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 11-265 amino acids of human ELF5 |
| Target Name | E74 like ETS transcription factor 5 |
| Other Names | ESE2 |
| Accession No. | Swissprot:Q9UKW6Gene Accession:BC029743 |
| Uniprot | Q9UKW6 |
| GeneID | 2001; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

The protein encoded by this gene is a member of an epithelium-specific subclass of the Ets transcription factor family. In addition to its role in regulating the later stages of terminal differentiation of keratinocytes, it appears to regulate a number of epithelium-specific genes found in tissues containing glandular epithelium such as salivary gland and prostate. It has very low affinity to DNA due to its negative regulatory domain at the amino terminus. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only