

Recombinant Human HHIPL2

Catalog No: #GP13573



Package Size: #GP13573-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HHIPL2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 196-451 amino acids of human HHIPL2
Target Name	HHIP like 2
Other Names	KIAA1822L
Accession No.	Swissprot:Q6UWX4 Gene Accession:BC015201
Uniprot	Q6UWX4
GeneID	79802;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only