

Recombinant Human C11orf54

Catalog No: #GP13581



Package Size: #GP13581-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C11orf54
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-169 amino acids of human C11orf54
Target Name	chromosome 11 open reading frame 54
Other Names	PTD012; PTOD012
Accession No.	Swissprot:Q9H0W9Gene Accession:BC012298
Uniprot	Q9H0W9
GeneID	28970;
Storage	-20~-80°C, pH 7.6 PBS

Background

Exhibits ester hydrolase activity on the substrate p-nitrophenyl acetate.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only