

Recombinant Human TTLL5

Catalog No: #GP13584



Package Size: #GP13584-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TTLL5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 39-223 amino acids of human TTLL5
Target Name	tubulin tyrosine ligase like 5
Other Names	STAMP; CORD19; KIAA0998
Accession No.	Swissprot:Q6EMB2Gene Accession:BC002766
Uniprot	Q6EMB2
GeneID	23093;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the tubulin tyrosine ligase like protein family. This protein interacts with two glucocorticoid receptor coactivators, transcriptional intermediary factor 2 and steroid receptor coactivator 1. This protein may function as a coregulator of glucocorticoid receptor mediated gene induction and repression. This protein may also function as an alpha tubulin polyglutamylase.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only