

## ZFAND2B Antibody

Catalog No: #43206

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	ZFAND2B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ZFAND2B protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human ZFAND2B
Target Name	ZFAND2B
Other Names	AIRAPL
Accession No.	Swiss-Prot#: Q8WV99?Gene ID: 130617
Uniprot	Q8WV99
GeneID	130617;
Calculated MW	28kd
Concentration	5.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:20-1:100

## Images



Gel: 12%SDS-PAGE  
Lysate: 40 µg  
Lane: Mouse brain tissue  
Primary antibody: 1/200 dilution  
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution  
Exposure time: 1 minute

## Background

This gene encodes a protein containing AN1-type zinc-fingers and ubiquitin-interacting motifs. The encoded protein likely associates with the

proteasome to stimulate the degradation of toxic or misfolded proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

---

Note: This product is for in vitro research use only