

Recombinant Human PDCL

Catalog No: #GP13591



Package Size: #GP13591-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDCL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-301 amino acids of human PDCL
Target Name	phosducin like
Other Names	PhLP
Accession No.	Swissprot:Q13371Gene Accession:BC017227
Uniprot	Q13371
GeneID	5082;
Storage	-20~-80°C, pH 7.6 PBS

Background

Phosducin-like protein is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin, a phosphoprotein expressed in retina and pineal gland. Both phosducin-like protein and phosphoducin have been shown to regulate G-protein signaling by binding to the beta-gamma subunits of G proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only