

Recombinant Human VTA1

Catalog No: #GP13592



Package Size: #GP13592-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VTA1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-307 amino acids of human VTA1
Target Name	vesicle trafficking 1
Other Names	DRG1; LIP5; SBP1; DRG-1; My012; C6orf55; HSPC228
Accession No.	Swissprot:Q9NP79Gene Accession:BC005937
Uniprot	Q9NP79
GeneID	51534;
Storage	-20~-80°C, pH 7.6 PBS

Background

C6ORF55 encodes a protein involved in trafficking of the multivesicular body, an endosomal compartment involved in sorting membrane proteins for degradation in lysosomes (Ward et al., 2005 [PubMed 15644320]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only