

## Recombinant Human SAP30BP

Catalog No: #GP13594



Package Size: #GP13594-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SAP30BP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-308 amino acids of human SAP30BP
Target Name	SAP30 binding protein
Other Names	HTRG; HTRP; HCNGP
Accession No.	Swissprot:Q9UHR5Gene Accession:BC007592
Uniprot	Q9UHR5
GeneID	29115;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only