

Recombinant Human PRKRA

Catalog No: #GP13599



Package Size: #GP13599-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PRKRA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-313 amino acids of human PRKRA
Target Name	protein activator of interferon induced protein kinase EIF2AK2
Other Names	RAX; PACT; DYT16; HSD14
Accession No.	Swissprot:O75569Gene Accession:BC009470
Uniprot	O75569
GeneID	8575;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only