

Recombinant Human SGTA

Catalog No: #GP13600



Package Size: #GP13600-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SGTA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-313 amino acids of human SGTA
Target Name	small glutamine rich tetratricopeptide repeat containing alpha
Other Names	SGT; hSGT; alphaSGT
Accession No.	Swissprot:O43765Gene Accession:BC000390
Uniprot	O43765
GeneID	6449;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only