

Recombinant Human DHRS7B

Catalog No: #GP13608



Package Size: #GP13608-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|-----------------------------------------------------------------------------------------|
| Product Name | Recombinant Human DHRS7B |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 1-325 amino acids of human DHRS7B |
| Target Name | dehydrogenase/reductase 7B |
| Other Names | CGI-93; SDR32C1 |
| Accession No. | Swissprot:Q6IAN0Gene Accession:BC004126 |
| Uniprot | Q6IAN0 |
| GeneID | 25979; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only