

Recombinant Human MSI2

Catalog No: #GP13609



Package Size: #GP13609-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human MSI2 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 129-328 amino acids of human MSI2 |
| Target Name | musashi RNA binding protein 2 |
| Other Names | MSI2H |
| Accession No. | Swissprot:Q96DH6Gene Accession:BC001526 |
| Uniprot | Q96DH6 |
| GeneID | 124540; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes an RNA-binding protein that is a member of the Musashi protein family. The encoded protein is transcriptional regulator that targets genes involved in development and cell cycle regulation. Mutations in this gene are associated with poor prognosis in certain types of cancers. This gene has also been shown to be rearranged in certain cancer cells.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only