

Recombinant Human TCF19

Catalog No: #GP13616



Package Size: #GP13616-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TCF19
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 146-345 amino acids of human TCF19
Target Name	transcription factor 19
Other Names	SC1; TCF-19
Accession No.	Swissprot:Q9Y242Gene Accession:BC002493
Uniprot	Q9Y242
GeneID	6941;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that contains a PHD-type zinc finger domain and likely functions as a transcription factor. The encoded protein plays a role proliferation and apoptosis of pancreatic beta cells. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only