

# Recombinant Human ENKD1

Catalog No: #GP13619



Package Size: #GP13619-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ENKD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 147-346 amino acids of human ENKD1
Target Name	enkurin domain containing 1
Other Names	C16orf48; DAKV6410
Accession No.	Swissprot:Q9H0I2Gene Accession:BC008284
Uniprot	Q9H0I2
GeneID	84080;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The exact function of ENKD1 remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only