

Recombinant Human CEP89

Catalog No: #GP13626



Package Size: #GP13626-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CEP89
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 144-343 amino acids of human CEP89
Target Name	centrosomal protein 89
Other Names	CEP123; CCDC123
Accession No.	Swissprot:Q96ST8Gene Accession:BC032307
Uniprot	Q96ST8
GeneID	84902;
Storage	-20~-80°C, pH 7.6 PBS

Background

Required for ciliogenesis. Also plays a role in mitochondrial metabolism where it may modulate complex IV activity.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only