

# Recombinant Human PPM1K

Catalog No: #GP13629



Package Size: #GP13629-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PPM1K
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 173-372 amino acids of human PPM1K
Target Name	protein phosphatase, Mg <sup>2+</sup> /Mn <sup>2+</sup> dependent 1K
Other Names	BDP; PTMP; PP2Cm; MSUDMV; PP2Ckappa; UG0882E07
Accession No.	Swissprot:Q8N3J5Gene Accession:BC037552
Uniprot	Q8N3J5
GeneID	152926;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the PPM family of Mn<sup>2+</sup>/Mg<sup>2+</sup>-dependent protein phosphatases. The encoded protein, essential for cell survival and development, is targeted to the mitochondria where it plays a key role in regulation of the mitochondrial permeability transition pore.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only