

# Recombinant Human SLC30A3

Catalog No: #GP13637



Package Size: #GP13637-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SLC30A3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 167-388 amino acids of human SLC30A3
Target Name	solute carrier family 30 member 3
Other Names	ZNT3
Accession No.	Swissprot:Q99726Gene Accession:BC028358
Uniprot	Q99726
GeneID	7781;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Involved in accumulation of zinc in synaptic vesicles.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only