

# Recombinant Human SIK3

Catalog No: #GP13638



Package Size: #GP13638-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SIK3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 971-1170 amino acids of human SIK3
Target Name	SIK family kinase 3
Other Names	L19; QSK; SEMDK; SIK-3
Accession No.	Swissprot:Q9Y2K2Gene Accession:BC035583
Uniprot	Q9Y2K2
GeneID	23387;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Inhibits CREB activity by phosphorylating and inhibiting activity of TORCs, the CREB-specific coactivators, like CRTC2/TORC2 and CRTC3/TORC3 in response to cAMP signaling.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only