

Recombinant Human DIAPH2

Catalog No: #GP13643



Package Size: #GP13643-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DIAPH2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 881-1080 amino acids of human DIAPH2
Target Name	diaphanous related formin 2
Other Names	DIA; POF; DIA2; DRF2; POF2; POF2A
Accession No.	Swissprot:O60879Gene Accession:BC117414
Uniprot	O60879
GeneID	1730;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene belongs to the diaphanous subfamily of the formin homology family of proteins. This gene may play a role in the development and normal function of the ovaries. Defects in this gene have been linked to premature ovarian failure 2. Alternatively spliced transcript variants encoding different isoforms have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only