

CASZ1 Antibody

Catalog No: #43214

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CASZ1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CASZ1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human CASZ1
Target Name	CASZ1
Other Names	CST; SRG; CAS11; ZNF693; dJ734G22.1
Accession No.	Swiss-Prot#: Q86V15 Gene ID: 54897
Uniprot	Q86V15
GeneID	54897;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43214 at dilution 1/20,

Background

The protein encoded by this gene is a zinc finger transcription factor. The encoded protein may function as a tumor suppressor, and single nucleotide polymorphisms in this gene are associated with blood pressure variation. Alternative splicing results in multiple transcript variants that encode different protein isoforms.

Note: This product is for in vitro research use only